***PLEASE***

***NOTE:***

*This document to be exported for proof and printing as pairs because the actual page size is that of two single pages i.e. 124 x 96mm*

# INSTRUCTIONS FOR USE <sup>1</sup>

Actifuse contains phase-pure, porous silicate substituted calcium phosphate granules for use as a synthetic bone graft substitute.

**Sterile - Single use - Store in a cool, dry place.**

## **Description:**

- 0.8 wt% silicate substituted calcium phosphate
- >99% phase pure silicate substituted calcium phosphate
- Total porosity 80%  $\pm$  2.5%
- Granule size 2 – 5 mm

**Indications:** as per morsellised grafting.

## **Directions for use:**

- Use in place of or added to morsellised allograft bone.
- Use in place of or added to morsellised autograft bone.
- Actifuse may be used alone to fill small unloaded bone defects, e.g. small tumours.
- Actifuse may be used alone or mixed with autograft or allograft bone for spinal fusion.
- Actifuse used alone may be mixed with either host blood or bone marrow.
- Mixing containers and stirrers must be sterilised before use.
- Actifuse implanted alone or mixed is used in the same manner as morsellised allograft.
- Actifuse may be used to fill an inter-body fusion cage.

**Risks and warnings:**

- The risks with Actifuse (e.g. infection, graft loosening, non-union, delayed union, repeat surgery) are the same as those for allograft bone. However, with Actifuse there is no risk of disease transmission.
- Actifuse should not be used where it could be subject to tension, torsion, compression, shear or bending. A conventional implant (e.g. screw, rod) can protect the graft from such load actions.
- Actifuse should not be used in volumetrically unconstrained sites (so the graft material cannot move or escape).
- Do not overfill the bony defect

**Sterility:**

- Sterilised by gamma irradiation.
- Do not re-sterilise.
- If the packaging or seal are damaged do not use. Discard.
- Discard unused product. Re-use of single-use devices creates a potential risk of patient or user infections.

**Contraindications:**

- Direct loading of graft material
- Impaction grafting for failed total hip or knee arthroplasty
- Patient metabolism may compromise bony regeneration
- Infection
- Inability to cover or deliberate non-coverage of graft site using soft tissue
- Avascular or compromised vascular network sites
- Excessive patient weight
- Medication which could slow bone healing

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